


PHOSPHORAMIDITES




Summit Pharmaceuticals Europe

Overview


Kolon has been supplying **refined high-purity APIs and Intermediates** using our core technology, **Sugar Chemistry**. Our sugar chemistry consists of (1) control of functional groups (2) various refinement & crystallization methods (3) scale-up experience

KLS Sugar / Nucleoside Chemistry




New Drug Intermediate / API


1. 2'-Modification derivatives



2. Base modification




3. Phosphorylation




Product List

Kolon has completed the synthesis of **23 Phosphoramidites** through sugar / nucleoside chemistry. Each of them can incorporate flow chemistry, scale up, specific quality control and customization.


DNA Phosphoramidites




TBDMS-RNA Phosphoramidites




2'-Modified Phosphoramidites (2'-OMe)



2'-Modified Phosphoramidites (2'-O-MOE)



2'-Modified Phosphoramidites (2'-F)



Quality Control

As the basic quality control methods for Phosphoramidites, Kolon conducts ^1H NMR, ^{31}P NMR, UV, TLC, HPLC, LC/MS, GC/MS and **MASS Spectrum** are used for stricter management of impurities that has adverse effects on Oligo synthesis.



^1H NMR

^{31}P NMR

UV –The UV test provides 4 values of data

TLC

HPLC
>98%

LC/MS

GC/MS

ESI-MS

Differentiation

- Securing product competitiveness through differentiated process development
- Preoccupation of the market through early establishment of GMP process

Stage 1. Fast Follower


- ✓ **Completion of general-purpose phosphoramidite development through competitor's product analysis**
 - Continuously expanding the product portfolio with customer requests
- ✓ **Promotion through collaboration with top 3 oligonucleotide makers around the world**

Stage 2. Differentiation

- ✓ **Securing cost competitiveness through differentiated process development**
 - Reduce manufacturing costs by developing switchable continuous manufacturing process suitable for small-scale multi-variety production
 - Development of crystallization process to replace column process
- ✓ **Diversification of Portfolio and development of customization**
 - Development of a product line that can be used to manufacture oligonucleotides in the future (Diester, Triester, etc.)
 - Development of various product lines to meet customer needs (Back-bone modification, 5'-capping, etc.)

Stage 3. First Mover

- ✓ **It is scheduled to build a production facility of 300 kg per year in the Kilo-Lab Center in 2025**
- ✓ **A production base of 3 MT per year will be established in the GMP plant in 2027**
 - There is a high probability that regulation will be strengthened as the RNA market grows



 **SPE**

Summit Pharmaceuticals Europe