


REAGENTS MADE IN JAPAN



 **SPE**

Summit Pharmaceuticals Europe



Hamari Chemicals, Ltd.

Unnatural Aminoacids

Extremely pure unnatural amino acids and peptides using Hamari's patented technology.



HODOGAYA CHEMICAL

1,1'-Carbonyldiimidazole

With exceptional stability and purity.



NIPPON CHEMICAL INDUSTRIAL

Buchwald Ligands

Highly effective phosphine ligands for Pd-catalyzed cross-coupling reactions.

(R,R)-QuinoxP*

(S,S)-QuinoxP*

P-chiral phosphine ligands with excellent enantioselectivity in asymmetric hydrogenations and C-C or C-N bond formation.



MANAC Incorporated

DAIB (Diacetoxyiodo)benzene

DAIM (Diacetoxyiodo)mesitylene

Low-toxic versatile oxidizing agents based on hypervalent iodine.

TOYOBO TOYOBO MC Corporation

Di-2-methoxyethyl azodicarboxylate (DMEAD)

A Mitsunobu reagent with easy work up.

Diphenylphosphoryl azide

Trimethylsilyl azide



S P E

Summit Pharmaceuticals Europe